

BTE SP & BTE UP: Oticon Xceed

| | Oticon Xceed 1 | Oticon Xceed 2 | Oticon Xceed 3 |
|---|--|--|--|
| Speech Understanding | | | |
| OpenSound Navigator™ | Level 1 | Level 2 | - |
| - Balancing power effect | 100% | 50% | - |
| - Max. noise removal | 9 dB | 5 dB | - |
| OpenSound Optimizer™ | • | • | • |
| Noise Reduction LX | - | - | • |
| Multiband Adaptive Directionality LX | - | - | • |
| Open Sound Booster | • | • | - |
| Speech Guard™ LX | Level 1 | Level 3 | - |
| Single Compression | - | - | • |
| Speech Rescue™ LX | • | • | • |
| Sound Quality | | | |
| Clear Dynamics | • | - | - |
| Spatial Noise Management | • | - | - |
| Processing Channels | 48 | 48 | 48 |
| Bass Boost (streaming) | • | • | • |
| Listening Comfort | | | |
| Transient Noise Management | 4 configurations | 3 configurations | - |
| Feedback Shield LX | • | • | • |
| Wind Noise Management | • | • | • |
| Personalisation & Optimising Fitting | | | |
| YouMatic™ LX | 3 configurations | 2 configurations | - |
| Fitting Bands | 14 | 12 | 8 |
| Oticon Firmware Updater | • | • | • |
| Multiple Directionality Options | • | • | • |
| Adaptation Management | • | • | • |
| VC Range and Step Size | • | • | • |
| Fitting Formulas | DSE/DSE linear, VAC+, NAL-NL1+2, DSL v5.0 | DSE/DSE linear, DSL v5.0, NAL-NL1+2, VAC+ | DSE/DSE linear, DSL v5.0, NAL-NL1+2, VAC+ |
| Connecting to the World | | | |
| Stereo streaming (2.4 GHz) | • | • | • |
| Made for iPhone® | • | • | • |
| Oticon ON App | • | • | • |
| ConnectClip, Remote Control 3.0, TV Adapter 3.0 | • | • | • |
| DAI/FM and EduMic Compatible | • | • | • |
| Special Requirements | | | |
| Bimodal Fitting Panel | • | • | • |
| CROS/BiCROS | • | • | • |
| Tinnitus SoundSupport™ | • | • | • |