

# miniRITE & miniBTE: Oticon Zircon™



	Oticon Zircon 1	Oticon Zircon 2
<b>Speech Understanding</b>		
OpenSound Navigator™	•	-
- Balancing power effect	40%	-
- Max. noise removal difficult/simple	6/0 dB	-
Multiband Adaptive Directionality	-	•
Noise Reduction	-	•
Speech Guard™	•	-
Single Compression	-	•
Frequency lowering	Speech Rescue™	Speech Rescue™
<b>Sound Quality</b>		
Fitting Bandwidth	8 kHz	8 kHz
Processing Channels	48	48
Bass Boost (streaming)	•	•
<b>Listening Comfort</b>		
Feedback Management	SuperShield & Feedback shield	SuperShield & Feedback shield
Wind Noise Management	•	•
Transient Noise Management	On/Off	-
<b>Personalisation &amp; Optimising Fitting</b>		
Fitting Bands	14	12
Multiple Directionality Options	•	•
Adaptation Management	•	•
Oticon Firmware Updater	•	•
Fitting Formulas	NAL-NL1+2, DSL v5.0	NAL-NL1+2, DSL v5.0
<b>Connecting to the World</b>		
Hands-free communication for iPhone	•	•
Direct streaming	•	•
Oticon ON & Oticon RemoteCare	•	•
ConnectClip, EduMic, Remote Control 3.0, TV Adapter 3.0, Phone Adapter 2.0	•	•
<b>Special Requirements</b>		
Tinnitus SoundSupport™	•	•
CROS/BiCROS Support	•	•